



Compact[®]
Peel



Compact Microdermabrasion System



Compact Peel Microdermabrasion

Over time, dead cells and cellular debris accumulate on the surface of the skin. This leads to a dull complexion and hinders the full absorption of cosmetics, thereby reducing the effectiveness of skin care products. Moreover, age produces fine lines, wrinkles and pigmentation spots that change our appearance. More effective than any other equipment of its kind on the market, the Compact Peel Microdermabrasion promotes a fresh radiant complexion and easily reduces the appearance of such imperfections.



Discover the Benefits

1 Innovative Treatment

The Compact Peel regenerating beauty cure is administered in a series of periodic treatments. Each treatment is performed in two phases: aspiration and microdermabrasion.

- **The aspiration phase** stimulates blood circulation to the skin's surface and supplies nutrients and energy to the epidermal cells.
- **The microdermabrasion phase** accelerates the regeneration process by promoting brand new, healthy and younger-looking skin. This skin shows greater consistency, helping to visibly reduce all signs of imperfections.

2 CristalVac Anti-Residue System

This unique procedure effectively removes the thin film of corundum powder left on the skin during microdermabrasion, resulting in a triple effect of deep exfoliation, regeneration and oxygenation of the epidermis. Thanks to CristalVac, your treatments will be safe, sanitary and effective.

3 DryFlo Anti-Moisture System

Very often, conventional microdermabrasion devices are prone to frequent blockage due to the build-up of condensation in the tubes, the pump and the applicator nozzle. Absorbent silica crystals allow the DryFlo Anti-Moisture System to remove all moisture circulating in the aspiration and propulsion systems, thereby preventing both the obstruction of microcrystals and the spread of bacteria during treatments.

4 Titanium Applicator Nozzle and Disposable Nozzle Caps

The Compact Peel's resistant titanium-alloy applicator nozzle withstands the wear of repeated use — even after several treatments. This allows it to retain its shape and consistently distribute an even umbrella jet spray. You benefit from the guarantee of quality treatment, without undue harshness or loss of crystals. The use of disposable, sealed and sterilized nozzle caps eliminates the risk of contamination and guarantees safe microdermabrasion treatments.

5 Corundum Microcrystals

Unlike chemical peels, microdermabrasion with corundum microcrystals is non-irritating and natural. Compact Peel uses corundum, a natural mineral second in hardness only to the diamond. The use of these extremely hard microcrystals is suitable for all skin types and allows you to offer exceptional microdermabrasion treatments. Corundum is also the most sterile and hypoallergenic product on the market.



Why Choose Compact Peel?

The Compact Peel is compact, effective and accessible, which makes it an indispensable device for every professional esthetician. It ensures that microdermabrasion treatments will give your clients youthful radiant skin — sure to be noticed. The Compact Peel consolidates two devices in one, and thereby allows two distinct treatment options. Additionally, the aspiration phase can enhance the effectiveness of your other beauty care treatments. The Compact Peel, which is easy to install and use, will allow you to perform effective microdermabrasion treatments with complete ease.

Perfection at your fingertips!

Features

Two distinct phases: Aspiration (Vacuum) and Microdermabrasion

DryFlo anti-moisture system to prevent blockage

CristalVac anti-residue system to remove crystals remaining on the skin

Motor with double pump

Rotary control for pressure adjustment

Wear-resistant titanium hand piece

Sterile disposable nozzle caps

Corundum microcrystals: 120 microns, extra fine molecule and hypoallergenic material

Pressure valve to clean all tubes

Device includes a 100-240 volt universal power supply

Meets international standards

Two-year warranty, training and advertising material available

Effects

Activates blood and lymph circulation

Stimulates cellular renewal

Eliminates dead skin cells and hyperkeratinization

Improves skin tone

Reduces wrinkles, fine lines and signs of ageing

Fades pigment spots and evens the complexion

Reduces cutaneous imperfections

Results

Before



After



Your Compact Peel Distributor:

Compact Peel is made by
Dectro International
1 800 463.5566
info@dectro.com
www.dectro.com



66-336115
090521

Printed in Canada