



- Thermolysis, thermocoagulation, Flash™ with 1/100 sec. precision.
- Improved MicroFlash™ power with 1/1000 sec. precision.
- OmniBlend™, Pulsing Blend™, electrolysis.
- "Processor Sentinel" protection circuit.
- Current stability monitor.
- MicroFlash™ Proportional Pulsing.
- "UL" key to calculate units of lye.
- HF energy level key "eL".
- Tolerance tester in Blend mode.
- 360-program guide for body area and hair coarseness (for non-insulated probes).
- Storage of 40 user-defined programs.
- Automatic detection for use without the footswitch.
- 9 cataphoresis and anaphoresis programs with alarm for each body area.
- Skin moisture tester.
- Treatment management software with cumulative income statements.
- 2 probeholder outlets.
- Optional external printer.
- Accessory tester.
- Optional iMM-Pac™ software.

Your Apilus® Distributor



apilus@dectro.com

Apilus®, the Gentle Way to Permanent Hair Removal™

Manufactured by Dectro International
www.dectro.com | 1.800.463.5566

Printed in Canada
070222



66-0123



APILUS^{MD}
SENIOR^{MC}TM

The Versatile Epilator



Give Some iMM-Pac™ to Your Treatments.



A High-Performance Epilator

The Apilus® Senior II™ boasts a 25% power increase in compared to its predecessor. This power, combined with adjustments accurate to one thousandth of a second ensure even more comfortable and rapid treatments with such innovations as the MicroFlash™ thermolysis mode. As well, the Apilus® Senior II™ offers a whole range of high-performance possibilities in electrolysis and Blend modes, combining different advantages to provide more efficient treatments. Take, for example, the unique OmniBlend™ mode, which can be combined with the Tolerance Test to ensure that treatments are perfectly adapted to your client's sensitivity.

Easy to Use

The Apilus® Senior II™ includes 360 treatment programs, carefully prepared for each specific area of the body, and conveniently displayed directly on screen. Furthermore, you can add up to 40 programs of your own: nothing could be simpler to adapt to each particular situation! As well, the addition of cataphoresis programs, combined with an application signal, greatly facilitates the use of this mode and allows you to end your epilation treatments easily and gently.

Innovative

Innovative techniques undeniably improve the quality of care provided by users of the Apilus® Senior II™. The MicroFlash™ Proportional Pulsing function automatically compensates the application time of thermolysis when more than one pulse is required, thus optimizing the comfort of your client. The Energy Level key (εL) indicates the calorific value applied during treatment in Flash™ or MicroFlash™ modes, confirming user-controlled adjustments. Furthermore, the Apilus® Senior II™ offers:

- An alphanumerical, multilingual display
- A smart adjustment of the insertion delay in automatic mode (without footswitch)
- An accessory test function
- A skin moisture analyser
- A simple billing system which includes the date, time, and the name and phone number of your salon
- Printed reports using the optional external printer

Performance ○

Versatility ○

Reliability ○

Simplicity ○

Innovation ○

Profitability ○

Profitable

Never has a management system been so easy to use. You simply select the option "begin of treatment", indicate the length of the session, and the Apilus® Senior II™ will do the rest! It records all the pertinent information related to your client's treatment, and displays the results on screen or prints them on the optional printer. You wish to compile your fees? A simple touch of a key will produce a printed profitability report. A useful and profitable system for your salon.

Adaptable

With the **iMM-Pac™** technology (optional), it is possible to add modes and programs to the ones that initially come with the Apilus® Senior II™. **iMM-Pac™** consists of a complementary application program which includes the MultiPlex™ and MultiBlend™ techniques and programs, the MicroFlash™ Proportional Pulsing – Variable Delay function, and a series of programs for insulated probes, for a total of 1080 programs. **Give some iMM-Pac™ to your treatments!**

Reliable

Taking advantage of the latest technological advances, the Apilus® Senior II™ represents a new level of excellence in computer-controlled permanent hair removal: comparison circuits guarantee unparalleled current stability, while the new "Processor Sentinel" protection module constantly monitors the condition of the microprocessor. It comes as no surprise that the Apilus® Senior II™ is the number one choice of professionals who insist on reliability, comfort and safety.

