

What is Lymphatic Circulation?

When thinking circulation, we very often focus on blood forgetting there is another important component flowing through the body: **lymph**.

It serves 3 purposes:

- 1) Evacuating excess fluids and toxins from tissues.
- 2) Absorbing and transporting fatty acids.
- 3) Producing immune-defense cells.

As opposed to the blood which is pumped by the heart, lymph only moves through muscular contractions. And it must circulate against the pull of gravity which doesn't make it any easier.

Without physical exertion, lymph stops moving and stagnates, promoting the accumulation of toxins and fluids in fat cells. When fluids and toxins become captive in these cells, trouble starts: you feel heaviness in the legs, cellulite appears and fatty deposits form on the hips and thighs.



Your Legs Deserve it!

Vita•Lymph™ Effects:


- Activates natural lymphatic and blood circulation
- Promotes the elimination of toxins and cellulite
- Prevents the appearance of varicose veins
- Soothes heaviness in the legs and edema
- Oxygenates tissues
- Accelerates metabolic reactions
- Gives a sensation of lightness and well-being

Vita•Lymph™ Technical Features:

- 16 independent pressure chambers working in pairs
- 9 customized programs for every type of circulation problem
- Special multiple program alternating between 4 different treatments
- Pressure: 1 to 20 kPa
- Even pressure in each chamber
- Ability to activate more than one chamber at a time
- Adjustable settings for a better-adapted treatment: treatment duration, pressure, and the time the pressure is maintained
- Can treat 4 areas simultaneously: feet, legs, groin, and abdomen
- Safety program preventing excess pressure in the pneumatic system
- Automatic pre-inflation at the beginning of the treatment to save time and improve results
- Automatic deflation at the end of the cycle
- Pause and Stop button to interrupt treatment at any time
- Inflatable legging system in 4 sections to adapt to any morphology
- Easy to use, handle and maintain
- No constant monitoring required
- Advanced design
- Large backlit screen

Your Vita•Lymph™ distributor

"Give your customer a moment of ultimate relaxation with this silent machine that will help them forget their day-to-day cares, improve their circulation and work at reducing the appearance of their cellulite."



Vita•Lymph
PRESSOTHERAPY

Made by Dectro International
1 800 463.5566
www.dectro.com



66-31621

080526
Printed in Canada

Vita•Lymph
PRESSOTHERAPY

Your Legs Deserve it!



Vita•Lymph™

To counter the negative effects of lymphatic circulation deficiency, *Vita•Lymph* pressotherapy treatments are just the thing because they help restore lymphatic and blood circulation, and eliminate excess fluids and toxins. The pressure exerted by *Vita•Lymph* during treatment simulates natural muscle movements which move the lymph in a manner similar to normal circulation.

Restoring lymphatic circulation:

- Reduces the feeling of heaviness in the legs.
- Decreases the appearance of cellulite, edema and fat deposits.
- Stimulates the immune-defense process by increasing production of white blood cells called lymphocytes.

Since an ounce of prevention is better than a pound of cure, *Vita•Lymph* also helps prevent the appearance of troubles associated with bad circulation. Your customers should not wait to suffer from all these problems before turning to *Vita•Lymph* pressotherapy. Stay ahead of their needs by offering preventative care they will be thankful for.

In addition to all these benefits, *Vita•Lymph* gives a feeling of lightness and well-being. After leaving the treatment chair, your customers will have the delicious feeling of floating on a cloud, the ultimate relaxation.

"Vita•Lymph gives a feeling of lightness and well-being."

Effective and Personalized

Vita•Lymph allows you to offer ultra effective pressotherapy care adapted to every customer. Thanks to its 16 air outlets working in pairs, it exerts pressure in more than one location at a time, from the abdomen to the soles of the feet, which manual drainage cannot provide. Sessions are relaxing, and the pressure felt gives a sensation similar to a massage.

To effectively care for every customer, you can choose one of 9 preset programs, each one targeting a specific imperfection. For even more customized care, you can create your own program by varying pressure values, treatment duration, and the time the pressure is maintained.

Also, *Vita•Lymph* allows you to combine up to 4 programs in the same session, which optimizes the treatment and improves its efficiency. You can, for instance, start with opening the lymph nodes, continue with a downward push from the abdomen to the soles of the feet, and finish with an upward push.

Advantageous

It is common for a professional to feel joint pain after doing a lot of manual drainage because the hands must work very vigorously. By using *Vita•Lymph* pressotherapy, there is no need to do manual massages, which makes the esthetician's work much easier. Additionally, since the pressure is exerted simultaneously on both legs, the sessions go faster and you can see more customers in a day.

Profitable

With the *Vita•Lymph* preset modes and safety system, you don't have to devote yourself entirely to the treatment, you can go do something else while the customer enjoys a moment of relaxation, or you can provide a complementary service such as a manicure or a facial. Your customers with busy schedules will appreciate the ability to get two treatments simultaneously. There's an easy way to make your *Vita•Lymph* profitable and increase the loyalty of your customers.

High-Performance!

Vita•Lymph pressotherapy care is the ideal choice since it applies to a wide customer base. It is even more interesting since it offers all the advantages of manual lymphatic drainage and more. *Vita•Lymph* technology is precise, easy to use, offers many possibilities, requires no effort or constant monitoring, and requires less application time.

In essence, this is a high-performance system!

