




What Is Cellulite?


Bad lymphatic circulation reduces the ability to eliminate toxins produced by the body. They then accumulate in fat cells.

When toxins accumulate in fat cells, a membrane of collagen fibers forms and sticks around these cells, preventing water and waste elimination. Unable to expel these elements, they inflate and agglomerate to form irregularly-shaped bunches under the skin surface producing the orange peel appearance.

Once cellulitic tissue is formed, the imperfection remains and resists even with exercise and diets. The most effective solution is to break the membrane and melt away the hardened connective tissue with the *Vital•Ultra* radio frequency which unblocks the subcutaneous fat. Additionally, using *Vital•Ultra* vacuodermie at the beginning and the end of the treatment facilitates the elimination of excess fluids and toxins by stimulating lymphatic circulation.



Cross-section of a skin without cellulite



Cross-section of a skin with cellulite

Perfect Silhouette

Vital-Ultra™ Effects

- Breaks down fat cells
- Promotes toxin elimination
- Reduces cellulite and fatty lumps
- Redefines contours
- Slims the silhouette (visible after the first treatment)
- Provides blood and lymphatic drainage
- Activates circulation and tissue oxygenation
- Eliminates swelling and bloating
- Accelerates metabolic processes
- Body remodeling (vacuodermie)
- Feeling of well-being and lightness
- Visible results right from the first treatment

Vital-Ultra™ Technical Features

- Leading-edge silhouette system
- Unique combination of radio frequency, infrared, light therapy, and vacuodermie producing impressive results
- 2 ergonomic handsets for radio frequency with infrared-light therapy
- 4 ultra-efficient radio frequency transmitters (3MHz)
- 16 electroluminescent diodes in the infrared spectrum (910nm)
- 8 electroluminescent diodes providing 4 colors: red, blue, green, yellow
- 5 treatment modes adapted to each customer: 1 lymphatic drainage mode (vacuodermie) + 4 cellulite modes (soft, hard, edematous, soft with capillaries)
- Massage with vacuum and roller for body remodeling (vacuodermie)
- For optimal care, you can add optional dermaphasic devices:
 - VitaPhase SM™ (StimuloMassage)
 - VitaPhase LT™ (21-diode light therapy, red, blue, green)
 - VitaPhase IR™ (infrared light therapy)
- Sophisticated and sleek design
- Large backlit screen
- Easy to use



Vital•Ultra
RADIO FREQUENCY - VACUODERMIE

Made by Dectro International
1 800 463.5566
www.dectro.com




080526
Printed Canada

Your *Vital•Ultra*™ distributor

Perfect Silhouette



Vital•Ultra
RADIO FREQUENCY - VACUODERMIE



**Vital•Ultra™,
the Solution for More Than 80% of Women!**

Many aspects of our lifestyle, such as stress, lack of physical activity, bad eating habits, and even high-heeled shoes, can impact blood and lymphatic circulation. Circulation problems are the main cause of cellulite, a problem affecting more than 80% of women. They are actually willing to pay considerable amounts to rid themselves of this horrible "orange-peel" skin.

Vital•Ultra represents the ultimate technological advance to visibly reduce the appearance of cellulite, redefine contours, slim the silhouette and improve circulation. It combines 4 technologies that are already established: radio frequency, infrared, light therapy, and vacuodermie.

In the field of esthetics, only the mixed technology of Vital•Ultra can effectively combat cellulite.

Radio Frequency


Each of the Vital•Ultra handsets has two 3-MHz radio frequency transmitters. Radio frequency is the perfect method to get rid of lumps of fat cells. The waves produced by Vital•Ultra create vibrations in the tissues that then turn into heat. Through the action of the vibration and heat, the membranes surrounding fat cells break apart and the cells are compressed, ejecting the captive toxins and water. These two elements are then absorbed by the lymphatic system whose action is also stimulated by heat.

After a Vital•Ultra radio frequency treatment, the appearance of cellulite is considerably reduced and the silhouette is slimmer.

Light Therapy and Infrared

The handsets of the Vital•Ultra also have 16 infrared diodes and 8 electroluminescent diodes emitting variable colors of red, yellow, green, and blue. The colored light stimulates the cells and the organism to correct imperfections. Simultaneously, infrared light dilates blood and lymphatic vessels, and regenerates elastic fibers in the dermis, which promotes circulation and reduces the risk of developing cellulite. The diodes, combined with radio frequency, enhance the effectiveness of the treatment.

Each color is associated to a specific treatment mode to perfectly target every type of cellulite.



Vacuodermie

An undeniable advantage of Vital•Ultra is that it offers vacuodermie. This technique performs a mechanical massage combining vacuum and the effects of a rolling manual massage. These movements emulate the effects of natural muscular contractions that allow lymphatic circulation in addition to transporting toxins released by radio frequency to the lymph nodes responsible for eliminating waste. Also, it causes the redistribution of subcutaneous fat which promotes the remodeling of contours.


Applied on its own or combined with radio frequency treatments, Vital•Ultra vacuodermie provides ultra-efficient drainage and remodeling, and consequently plays a major role in eliminating and preventing cellulite.

Customized Treatments

Thanks to Vital•Ultra's specific programs, you can treat each customer according to their type of cellulite and therefore provide customized, optimal and progressive care.

Vital•Ultra™ Specific programs

Program 1	Hard Cellulite	Pulsed 3-MHz radio frequency + 910 nm infrared emission modulated to 5 KHz + red light	
Program 2	Edematous Cellulite	Continuous 3-MHz radio frequency + 910 nm infrared emission modulated to 3.5 KHz + yellow light	
Program 3	Soft Cellulite	Continuous 3-MHz radio frequency + 910 nm infrared emission modulated to 4.5 KHz + green light	
Program 4	Soft Cellulite with Capillaries	Rapidly pulsed 3-MHz radio frequency + 910 nm infrared emission modulated to 3.5 KHz + green and blue lights	



Before After

5 Vital•Ultra™ treatments

Vital•Ultra allows you to offer your customers a unique combination of anti-cellulite care through radio frequency, infrared, light therapy, and vacuodermie. These treatments contribute to reducing cellulite, redefining contours, slimming down the silhouette, eliminating swelling, and providing a feeling of wellbeing. Right from the first treatments, your customers will notice a reduction of a few centimeters in the treated areas.

Show your customers the wonders of the Vital•Ultra silhouette care!