

# *VitaPeel* <sup>ion</sup> MC TM

*Crystallize your beauty*



For effective and comfortable  
**microdermabrasion,**  
opt for the advanced VitaPeel™ Ion system!

The VitaPeel™ Ion system has revolutionized the microdermabrasion market with its Ion Mode and the VitaPhase™ IM and SM. Thanks to these functions, found only on the VitaPeel™ Ion, you now have the ability to offer complementary iontophoresis, stimulomassage and ionomassage treatments to your clientele!

### Ion Mode

Through the use of a galvanic current, the Ion mode allows you to perform specific iontophoresis treatments (ionization/deep-cleansing). Ionization enhances the absorption of ionizable cosmetic products, while the deep-cleansing function, thanks to its hyperemic and vasodilating effect, increases the skin's permeability and yields a truly deep cutaneous cleansing. These high-performance treatments can be performed using either the VitaPhase™ IM device or the multipurpose electrodes, depending on the client's specific needs.

### VitaPhase™ IM

The VitaPhase™ IM (ionomassage) allows you to perform iontophoresis treatments as well as active vibro-massages in synergy with the Ion mode. It offers two distinct modes: iontophoresis and combined iontophoresis/vibro-massage. After a microdermabrasion session, combining the VitaPhase™ IM with an ionizable soothing gel will help reduce redness, relieve the sensation of heat and restore the pH level. Furthermore, performed in conjunction with a vibro-massage, iontophoresis promotes an even greater penetration of active products, thereby increasing the effectiveness of your skincare treatments.

### VitaPhase™ SM

The VitaPhase™ SM (stimulomassage), on the other hand, is designed to perform cellular renewal treatments, facial stimulation with pulsing current as well as relaxing vibro-massages. Equipped with 10 stimulation electrodes and 5 light-emitting diodes, it offers two distinct modes: relaxing vibro-massage and combined facial stimulation/vibro-massage. In synergy with these two modes, the diodes' red light emission has a chromotherapeutic effect that revitalizes. After a microdermabrasion session with the VitaPeel™ Ion, the best course of action is to relieve tension with a relaxing vibro-massage and stimulate the cellular renewal process in order to obtain the best possible results.



For effective and comfortable  
**microdermabrasion,**  
opt for the advanced VitaPeel™ Ion system!

