






APILUS[®]
[COMPUTERIZED PERMANENT EPILATION]
[ÉPILATION DÉFINITIVE INFORMATISÉE]

C · L · É · O[™]



The Simple
Computerized Epilator

Hair Removal by Cléo, simple

You're starting out in computer-assisted permanent electroepilation and are searching for a simple yet upgradeable electrolysis device? Opt for the Apilus Cléo, it's simple and can evolve right along with you! Available in two versions, the Cléo allows you to move from a basic to a more advanced version at any time. A simple way to meet your current and future needs. Brilliant!



1

Start with the Flash 128 version

An easy-to-use electroepilation system that yields permanent results.

- Equipped with world-renowned, reliable and ultra stable electronic circuitry.
- Uses 13.56 MHz rapid frequency and digital technology for precise and effective treatments.
- Includes Thermolysis and Flash, two comfortable techniques that quickly eliminate all types of hair.
- Offers 128 user-friendly programs that guide you through the adjustment of your parameters, for all body areas and hair types.
- Includes an automatic function to convert programs for use with insulated probes, so you don't have to think about reducing your settings.

2

Move on to the Blend 256 version

A flexible system that offers even more options and treatment modes.

- Includes all Flash 128 version features and more!
- Adds the simplified Blend mode to its functions, the ideal solution for distorted hairs.
- Offers the Tolerance Test, an essential function to set your Blend treatments according to your clients' sensitivity.
- Adds a range of Blend programs, for a total of 256, to facilitate your work.
- Integrates the Cataphoresis post-epilatory mode, to soothe the skin after treatment.

1 system. 2 versions. Simple.



The Gentle Way to Permanent Hair Removal

Your Apilus® Distributor



66-1378

110413

Printed in Canada

Manufactured by Dectro International

1 800.463.5566 • apilus@dectro.com • www.dectro.com