Cross Reference Chart APILUS



		xCell Pro	xCell Pur	MILLA	SENIOR 3G	CENTA	JUNIOR 3G
Frequency	13.56 MHz			✓	✓	✓	✓
	27.12 MHz	✓	✓				
Power		·····	·····	····	···	··	·
Speed	Milliseconds	✓	✓	✓	✓		
	Centiseconds					✓	✓
Thermolysis	PicoFlash (ultrafast 27.12 MHz HF current)	✓	✓				
	MeloFlash (fast 27.12 MHz HF current)	✓	✓				
	Synchro (synchronized application of multiple ultrafast 27.12 MHz HF micropulses)	✓	optional				
	MultiPlex (slow thermolysis + MicroFlash or PicoFlash pulse at the end)	✓	✓	✓	✓		
	MicroFlash (very fast 13.56 MHz HF current)			✓	✓		
	Flash (fast 13.56 MHz HF current)			✓	✓	✓	✓
	Thermolysis (13.56 MHz HF current)					✓	✓
	Pp/vd function (allows the application of two reduced MicroFlash or PicoFlash pulses in two different spots inside the hair follicle)	✓	✓	✓	✓		
Blend	OmniBlend (GLV + HF currents)	✓	✓	✓	✓	✓	✓
	MultiBlend (GLV + HF currents + MicroFlash or PicoFlash pulse at the end)	✓	✓	✓	✓		
	SynchroBlend (continuous GLV current + multiple 27.12 MHz HF micropulses)	✓	optional				
	EvoluBlend (simultaneous application of a gradually increasing GLV current and 27.12 MHz HF current)	✓	optional				
	PicoBlend (continuous GLV current + 27.12 MHz HF pulses)	✓	✓				
	Pulsing Blend (continuous GLV current + 13.56 MHz HF pulses)			✓	✓	✓	✓
Electrolysis (continuous GLV current)		✓	✓	✓	✓	✓	✓
Printer		✓					
Tolerance test		✓	✓	✓	✓	✓	✓
Moisture test*		✓	✓		✓		✓
Number of programs (covering each body area, hair size and probe type)		29 700	19 800	18 744	1080	11 044	512
Programs based on the client's hormonal profile, type of isolated filaments, and filament size.		✓	✓	✓		✓	
Number of personalized programs		100	60		60		
Cataphoresis mode (soothing treatment)		✓	✓	✓	✓	✓	✓
Anaphoresis mode		✓	✓	✓	✓	✓	✓
VitaPhase Ionomassage* device		✓	optional		optional		optional
Screen	Regular (80 characters)				✓		✓
	10-inch touchscreen (25.5 cm)			✓			
	7-inch touchscreen (18 cm)	✓	✓			✓	
Key to calculate U.L. (units of lye)		✓	✓	✓	✓	✓	
Key to calculate εL (energy level)		✓	✓	✓	✓	✓	
Operation mode without footswitch		✓	✓	✓	✓	✓	✓
Management software and 1000 client files		✓					
Medical records, appointment calendar and statistics		✓					
Smart card system		✓					
USB port and updates		✓	✓	✓		✓	
Number of probe holder outlets		2	2	2	2	1	1
Accessory test		✓	✓	✓	✓		
Power test				✓		✓	
Current protection and stability monitor		✓	✓	✓	✓	✓	✓
International power supply		✓	✓	✓		✓	

* not available in the USA